



MODEL: TR2001 - GTE

Basic Features:



HI-LO (605-1010mm) elevation by HYDRAULIC FOOT PEDAL



Backrest (75°) elevation operated by GAS LIFT



Trendelenburg (12°) by MANUAL GAS ACTIVATION



Reverse Trendelenburg (6°) by MANUAL GAS ACTIVATION

Standard Accessories:

1. Stainless steel IV Drip Pole
2. 50mm mattress
3. Oxygen holder
4. 200mm Swivel castor with directional lock and central lock system

## TREATMENT / RECOVERY TROLLEY

### MODEL: TR2001 - GTE

|               |  |
|---------------|--|
| MATERIAL      | : MILD STEEL FRAMEWORK                           |
| TREATMENT     | : ZINC-PHOSPHATE WASHING                         |
| FINISHING     | : EPOXY POWDER COATING                           |
| DIMENSION     | : L2145 x W890 x H605-1010mm                     |
| NETT WEIGHT   | : 130kg  |
| LOAD CAPACITY | : < 250kg  |
| CASTOR DIA.   | : Ø8" / 200mm swivel castor with central locking |

#### FEATURES :

- Backrest (75°) elevation operated by gas lift
- Hi-Lo elevation operated by hydraulic foot pedal
- 2-way Tilting – Trendelenburg (12°) & Reverse Trend (6°) by gas activation
- Collapsible PU handle bar, at both end of trolley
- Collapsible Epoxy-finish safety side rails with latch lock (Dim : L1525 x H415mm)
- 50mm foam mattress with PVC covering
- Oxygen tank holder ('E' size)
- Stainless Steel IV Pole with double hook
- 4 IV pole attachment
- Ventilated 2-panel mattress base
- Roller bumpers at all corners, to protect body frame
- IV Pole (IVP 1)



Foot pedal to activate central lock, directional, free of wheels



Foot end lever to operate Trendelenburg & Reverse Trendelenburg

