



MODEL: MET 6600 I / MET 6600 II


Basic Features:




HI-LO (630-1030mm) elevation by HYDRAULIC FOOT PEDAL



Backrest (75°) elevation operated by GAS LIFT



Trendelenburg (18°) by FOOT PEDAL OPERATED



Reverse Trendelenburg (18°) by FOOT PEDAL OPERATED

Standard Accessories:

1. Stainless steel IV Drip Pole
2. 50mm mattress
3. Oxygen holder
4. 200mm Swivel castor with directional lock and central lock system

Optional Accessories:

1. Removable patient monitor / Defibrillator tray

MODEL: MET 6600 I / MET 6600 II

MATERIAL	: MILD STEEL FRAMEWORK
TREATMENT	: ZINC-PHOSPHATE WASHING
FINISHING	: EPOXY POWDER COATING
DIMENSION	: L2130 x W800 x H560~900mm
LOAD CAPACITY	: <300kg
CASTER DIAMETER	: 200mm (Central Locking)

FEATURES :

- Backrest adjustment by Gas Lift assisted (0~90°)
- Legrest elevation by manual hand crank (0~40°) (MODEL: MET 6600II)
- Foot adjustment by manual ratchet (MODEL: MET 6600II)
- Dual hydraulic Hi-Lo elevation operated by foot pedal on both sides
- Trendelenburg (18°) and Reverse Trendelenburg (18°) adjustment by handsfree foot operation on both sides
- 200mm swivel castors with 5th wheel
- Central Brake & Directional steer pedal at head, foot, right and left of trolley
- Collapsible ¾ length Chromed-finish safety side rails with spring-loaded lock (Dim : L1400 x H340mm)
- Foldable handles at head & foot end of bed
- ABS mattress base with restraint strap holders for easy cleaning purposes.
- 75mm thick foam mattress
- Seamless welded mattress cover
- Roller bumpers with protection cover on all corners of trolley
- Vertical oxygen tank holder ('E' size)
- 2 Mounted Stainless Steel IV Pole with double hook
- Utility Tray
- Drainage bag holder bars on both sides of bed



Backrest activated by gas lift & Caster Locking at 4 sides of bed



To operate HILO, Trendelenburg or Reverse Trendelenburg by foot pedal